


## **FodiCuff**<sup>TM</sup>

Reusable Blood Pressure Cuffs



## **FodiCuff™**

#### Reusable Blood Pressure Cuffs





#### Reusable Blood Pressure Cuffs

If you need a reusable blood pressure cuff, you have found what you are looking for. Our blood pressure cuffs, made of durable and long-lasting materials, appeal to a wide audience and have been developed to adapt to all difficult conditions that may occur in the clinical environment.


**Easy Application and Placement:** Thanks to its improved design, it facilitates application to the patient. Color-coded dimensions and range markings have been prepared considering all the details required for a perfect fit and accurate measurement.

**Durability:** Our blood pressure cuffs, made of water-repellent polyester fabric, are long-lasting and resistant to water etc. It ensures minimum cross-contamination as it removes the substances. We offer you our velcro fasteners, which are made of polyamide and polyester and have a minimum opening-closing life of 10,000 times, supported by professional stitches and first-class materials.

**100% Reliable Measurement:** Each reusable blood pressure cuff is delivered to you after passing the leakage test for accurate measurements and values.

**Quality Standards:** FODI Blood Pressure Cuff is developed and produced in accordance with the EU 2017/745 Medical Devices Regulation, in CE and TSE standards. Multiple Configurations: With various connector options, our cuffs are compatible with devices produced by all major healthcare facility equipment manufacturers in the world and also offer single and double hose options.

**Material:** For the safety of the patient, we do not use any materials containing latex, PVC or DEHP, and we produce our inner cuffs from high pressure resistant TPU material.



### **FodiCuff™**

#### Reusable Blood Pressure Cuffs





#### Size of FodiCuff™

Symbol	Size	Designation	Cuff	Bladder	Arm Circumference	Single Hose	Double Hose
ŕ	13	Thigh	19 x 83 cm	21 x 40 cm	46 x 66 cm	RBC - S06	RBC - D06
Ť	12	Large Adult	17 x 62 cm	16 x 35 cm	33 x 47 cm	RBC - S05	RBC - D05
Ť	11	Adult	18 x 47 cm	13 x 27 cm	25 x 35 cm	RBC - S04	RBC - D04
Ť	10	Small Adult	12 x 41 cm	10 x 21 cm	20 x 26 cm	RBC - S03	RBC - D03
Ť	9	Child	8 x 27 cm	7 x 15 cm	14 x 20 cm	RBC - S02	RBC - D02
å	8	Infant	6 x 21 cm	6 x 9 cm	10 x 14 cm	RBC - S01	RBC - D01
~	7	Small Infant	5 x 17 cm	5 x 7.5 cm	6 x 10 cm	RBC - S08	RBC - D08
<b>3</b> 6	6	Newborn	4.5 x 14 cm	4 x 6 cm	4 x 8 cm	RBC - S07	RBC - D07

#### **BP Cuff Connectors**





Female Quick Connect Double Hose / CS01-CS01



Male Quick Connect Double Hose / CL01-CL01





Male Quick Connect Single Hose / CL01



Female Quick Connect Single Hose / CS01



Threaded (Screw) Single Hose / MA01

#### **Symbols**



Symbol indicating Caution



Index Line



See Instruction for use



Not made with natural rubber latex



Product in compliance with (EU)2017/745 Medical Devices Regulation



Artery symbol and arrow should be placed over brachial or femoral artery

#### **Design Features**

Cuff Components	Material
Cuff Shell	Nylon,PES Waterproof Polyester
Hook	Nylon
Loop	Nylon
Theread	Poliamid
Bladder	TPU (Thermoplastic Polyurethane
Hose	TPE, TPU







# **FODI MEDICAL®**

#### **FOD Group Industrial Trade Limited Company**

TEL: +(90) 312 3941118 E-mail: fodi@fodimedical.com

export@fodgrup.com WEBSITE: www.fodimedical.com ADD: Ivedik Industrial Zone 1358. Street

No:30 06378

Yenimahalle, Ankara, Turkey



